

LABINCO L27 Large Rotary Mixer



Technical data

Number of mix positions	1 bar, drum, bottle holder
Mixing quantity max.	4 kg
Motor rating input	85 W
Motor rating output	55 W
Speed display	analog
Speed range	1 - 50 rpm
Timer display	digital 1s. up to 99h 59m 59s
Heat output	-
Heating rate	-
Heating temperature range	-
Heat control	-
Heat control accuracy	-
Speed control	stepless
Adjustable safety circuit min.	-
Adjustable safety circuit max.	-
Connection for ext. temperature sensor	-
Control accuracy with sensor	-
Temperature constancy in medium	-
Set-up plate material	Steel/nylon
Set-up spring clip dimensions	see www.labinco-bv.com
Dimensions (W x D x H)	540 x 160 x 235 mm
Weight	6,2 kg
Permissible ambient temperature	5 - 40 °C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
RS 232 interface	-
Product number	27000
Voltage	220 - 230 V (115 / 100 V)
Frequency	50/60 Hz
Power input	100 W

LABINCO L27



Large Rotary Mixer



Discription ; The Model LD27 "ROTARY MIXER WITH DIGITAL TIMER" is a instrument that provides Bottle accesory/Hybridisation accesory/Top over bottom accesory or Extraction accesory. Provided with a strong dc motor up to 0-50 rpm. Timer from 1 sec up to 99 h-59 min-59 sec.

Operating instructions ; Unpack the Model LD27 and choose a flat surface table. After choosing proper table-space, unpack assesory and place it between motorside-support and handgripside support. Put now the linecord into a suitable wall contact (see specifications). You can choose to continu or timer function. For using the continuous function press on START (or STOP) button, for using the timer see "Setting the time" below here.

Key Functions; (PROG) Selects between hours, minutes and seconds. (+) Increments the selected time range. (-) Decrements the selected time range. (RESET) Resets to start time. (START/STOP) Starts and stops the countdown.

Setting the time; Press the (PROG) button once. This will cause the hour digit to flash. Using the (+) (-) keys to set the desired amount of hours required. Repeat the process to set minutes and seconds. Pressing the (START/STOP) button activates the relay and the set time begins to count down in one second steps. When the set time reaches zero the unit emits an audible signal, the relay deenergizes and the set time is indicated on the display.

LABINCO BV

P.O.BOX 3259 4800DG BREDA THE NETHERLANDS

ZINKSTRAAT 64 4823AC BREDA THE NETHERLANDS

TEL;0031(0)76-5417582

FAX;0031(0)76-5420173

[E-mail;info@labinco-bv.com](mailto:info@labinco-bv.com)

<http://www.labinco-bv.com>