


LABINCO LD79 Digital Test Tube Rotator

New, from 1 rpm up to 99 rpm!!!



Description

The "Labinco model LD79" testtube rotator is an instrument that provides serial inclined isolation of testtubes, with diameter-range from 4 up to 45 mm. All powder coated and stainless steel parts. Provided with stepping motor 1-99 rpm, permites tube exchange while operating. Angle of inclined rotation can be continously set between 90 and 27 degr. Equipped with timer up to 99 hrs, 99 min. Comes with line-cord;but without discs. Standart model will operate on 230 volt, 50/60HZ. (other lines on request).

Operating instructions

Unpack Model LD79 and carefully remove all packing material. After choosing proper table space you can make proper line-voltage connection to (see specifications). You can start model ld 79 by switching on power switch and setting the speed, press start and the motor will run as you can see from the drive nut. Choose the disc you need from the list. Unscrew the black plastic nut from the drive spindle,place the disc in the proper fixed position and screw on the plastic nut. Tighten this nut rather strongly. You can choose the rotating angle by loosening the black tightening wheels on both sides and turning the motor-block to the right angle. Fasten the black tightening wheels. Now the unit is ready for use.

LABINCO BV

P.O.BOX 3259 4800DG BREDA THE NETHERLANDS

ZINKSTRAAT 64 4823AC BREDA THE NETHERLANDS

TEL;0031(0)76-5417582


FAX;0031(0)76-5420173

[E-mail;info@labinco-bv.com](mailto:info@labinco-bv.com)

<http://www.labinco-bv.com>

LABINCO LD79 Digital Test Tube Rotator

New, from 1 rpm up to 99 rpm!!!



Technical data

Number of mix positions	1 for max 2 disks
Mixing quantity max.	1 kg
Motor rating input	50 W
Motor rating output	35 W
Speed display	digital
Speed range	1 - 99 rpm
Timer display	digital 1s. up to 99h 59m 59s
Heat output	-
Heating rate	-
Heating temperature range	-
Heat control	-
Heat control accuracy	-
Speed control	digital
Adjustable safety circuit min.	-
Adjustable safety circuit max.	-
Connection for ext. temperature sensor	-
Control accuracy with sensor	-
Temperature constancy in medium	-
Set-up plate material	Steel/nylon clamps
Set-up spring clip dimensions	see www.labinco-bv.com
Dimensions (W x D x H)	200 x 200 x 250 mm
Weight	3.5 kg
Permissible ambient temperature	5 - 40 °C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
RS 232 interface	-
Product number	79000
Voltage	220 - 230 V (115 / 100 V)
Frequency	50/60 Hz
Power input	60 W