


LABINCO L51 12-channel fixed speed pump



Technical data

Number of tubes positions	12
Pump quantity max.	0.03-3.70 ml/min
Motor rating input	75W
Motor rating output	55 W
Speed display	-
Speed range	10 rpm not adjustable
Timer display	-
Heat output	-
Heating rate	-
Heating temperature range	-
Heat control	-
Heat control accuracy	-
Speed control	not adjustable
Adjustable safety circuit min.	-
Adjustable safety circuit max.	-
Connection for ext. temperature sensor	-
Control accuracy with sensor	-
Temperature constancy in medium	-
Set-up plate material	6 x Stainless steel rolls
Set-up dimensions	6 x rolls Ø 7.8 mm
Dimensions (W x D x H)	260 x 130 x 120 mm
Weight	7 kg
Permissible ambient temperature	5 - 60 °C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
RS 232 interface	-
Analog output	-
Voltage	220 - 230 V (115 / 100 V)
Frequency	50/60 Hz
Power input	75 W