


LABINCO LD209450 Rolling Bench

450 mm long rolls, Ø 46 mm.



Technical data

Number of rolling positions	1 up to 4
rolling quantity max.	10 kg
Motor rating input	75W
Motor rating output	50 W
Speed display	digital
Speed range	1-99 rpm
Timer display	digital 1s. up to 99h 59m 59s
Heat output	-
Heating rate	-
Heating temperature range	-
Heat control	-
Heat control accuracy	-
Speed control	digital
Adjustable safety circuit min.	-
Adjustable safety circuit max.	-
Connection for ext. temperature sensor	-
Control accuracy with sensor	-
Temperature constancy in medium	-
Set-up plate material	Steel with PE coated
Set-up dimensions	12 up to 2150 mm
Dimensions (W x D x H)	650 x 350 x 130 mm
Weight	9.4 kg
Permissible ambient temperature	5 - 50 °C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
RS 232 interface	-
Product number	20904
Voltage	220 - 230 V
Frequency	50/60 Hz
Power input	100W

